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Method		Other PCP Rapid Test		Total Results
PCP Test	Results	Positive	Negative	
	Positive	55	0	55
	Negative	2	155	157
Total Results		57	155	212
Agreement with this Rapid Test Kit		97%	100%	99%

When compared at 25 ng/mL cut-off with GC/MS, the following results were tabulated:

Method		GC/MS		Total Results
PCP Test	Results	Positive	Negative	
	Positive	50	5	55
	Negative	0	157	157
Total Results		50	162	212
Agreement with this Rapid Test Kit		100%	97%	98%

Analytical Sensitivity

A drug-free urine pool was spiked with Phencyclidine at the following concentrations: 0 ng/mL, 12.5 ng/mL, 18.75 ng/mL, 25 ng/mL, 31.25 ng/mL and 37.5 ng/mL. The result demonstrates >99% accuracy at 50% above and 50% below the cut-off concentration. The data are summarized below:

PCP Concentration (ng/mL)	Percent of Cut-off	n	Visual Result	
			Negative	Positive
0	0	30	30	0
12.5	-50%	30	30	0
18.75	-25%	30	26	4
25	Cut-off	30	11	19
31.25	+25%	30	8	22
37.5	+50%	30	0	30

Analytical Specificity

The following table lists the compound that is positively detected in urine by PCP Urine Screening test cassette at 5 minutes.

Compound	Concentration (ng/mL)
4-HydroxyPhencyclidine	12,500
Phencyclidine	25

Precision

A study was conducted at 3 physicians' offices by untrained operators using 3 different lots of product to demonstrate the within run, between run and between operator precision. An identical panel of coded specimens containing, according to GC/MS, no Phencyclidine, 25% Phencyclidine above and below the cut-off, and 50% Phencyclidine above and below the 25 ng/mL cut-off was provided to each site. The results are tabulated:

Phencyclidine conc. (ng/mL)	n per site	Site A		Site B		Site C	
		-	+	-	+	-	+
0	15	15	0	15	0	15	0
12.5	15	15	0	15	0	15	0
18.75	15	14	1	12	3	14	1
31.25	15	1	14	6	9	6	9
37.5	15	1	14	0	15	0	15

Effect of Urinary Specific Gravity

Fifteen (15) urine specimens with specific gravity ranging from 1.001 to 1.032 were spiked with 12.5 ng/mL and 37.5 ng/mL of Phencyclidine respectively. PCP Urine Screening test cassette was tested in duplicate using the fifteen neat and spiked urine specimens. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

Effect of Urinary pH

The pH of an aliquoted negative urine pool was adjusted to a pH range of 5 to 9 in 1 pH unit increments and spiked with Phencyclidine to 12.5 ng/mL and 37.5ng/mL. The spiked, pH-adjusted urine was tested with PCP Urine Screening test cassette in duplicate. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

Cross-Reactivity

A study was conducted to determine the cross-reactivity of the test with compounds in either drug-free urine or Phencyclidine positive urine. The following compounds show no cross-reactivity when tested with PCP Urine Screening test cassette at a concentration of 100 µg/mL.

Non Cross Reacting Compounds

Acetaminophen	β-Estradiol	Oxymetazoline
Acetophenetidin	Estrone-3-sulfate	Papaverine
N-Acetylprocainamide	Ethyl-p-aminobenzoate	Penicillin-G
Acetylsalicylic acid	Fenoprofen	Pentazocine hydrochloride
Aminopyrine	Furosemide	Pentobarbital
Amitypyline	Gentisic acid	Perphenazine
Amobarbital	Hemoglobin	Phenelzine
Amoxicillin	Hydralazine	Phenobarbital
Ampicillin	Hydrochlorothiazide	Phentermine
L-Ascorbic acid	Hydrocodone	L-Phenylephrine
D,L-Amphetamine	Hydrocortisone	β-Phenylethylamine
Apomorphine	O-Hydroxyhippuric acid	Phenylpropanolamine
Aspartame	p-Hydroxy-methamphetamine	Prednisolone
Atropine		Prednisone
Benzilic acid	3-Hydroxytyramine	Procaine
Benzoic acid	Ibuprofen	Promazine
Benzoylcegonine	Imipramine	Promethazine
Benzphetamine	Iproniazid	D,L-Propranolol
Bilirubin	(±) - Isoproterenol	D-Propoxyphene
(±) -Brompheniramine	Isoxsuprine	D-Pseudoephedrine
Caffeine	Ketamine	Quinidine
Cannabidiol	Ketoprofen	Quinine
Cannabinal	Labetalol	Ranitidine
Chloralhydrate	Loperamide	Salicylic acid
Chloramphenicol	Maprotiline	Secobarbital
Chlordiazepoxide	Meperidine	Serotonin (5-Hydroxytyramine)
Chlorothiazide	Meprobamate	
(±) Chlorpheniramine	Methadone	Sulfamethazine
Chlorpromazine	Methoxyphenamine	Sulindac
Chlorquine	(+)	Temazepam
	3,4-Methylenedioxy-amphetamine	
Cholesterol		Tetracycline
Clomipramine	(+)3,4-Methylenedioxy-	Tetrahydrocortisone,

Clonidine	methamphetamine	3 acetate
	Morphine-3-β-D glucuronide	Tetrahydrocortisone
	Morphine Sulfate	3 (β-D glucuronide)
Cocaine		Tetrahydrozoline
hydrochloride		
Codeine	Nalidixic acid	Thiamine
Cortisone	Naloxone	Thioridazine
(-)-Cotinine	Naltrexone	D, L-Tyrosine
Creatinine	Naproxen	Tolbutamide
Deoxycorticosterone	Niacinamide	Triamterene
Dextromethorphan	Nifedipine	Trifluoperazine
Diazepam	Norcodein	Trimethoprim
Diclofenac	Norethindrone	Trimipramine
Diflunisal	D-Norpropoxyphene	Tryptamine
Digoxin	Noscapine	D, L-Tryptophan
Diphenhydramine	D,L-Octopamine	Tyramine
Doxylamine	Oxalic acid	Uric acid
Ecgonine	Oxazepam	Verapamil
hydrochloride		
Ecgonine methylester	Oxolinic acid	Zomepirac
(-)-Y Ephedrine	Oxycodone	
Erythromycin		

BIBLIOGRAPHY

- Robert DeCresce. Drug Testing in the workplace. 114
- Baselt RC. Disposition of Toxic Drugs and Chemicals in Man. 2nd Ed. Biomedical Publ., Davis, CA. 1982; 488
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986