AB0 set MP

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The diagnostic set is used at the patient's bedside immediately before the transfusion for direct examination of compatibility between the recipient blood and the donor blood in the AB0 blood group system.

General information

The diagnostic set for AB0 blood grouping (monoclonal) contains one vial of monoclonal Anti-A reagent and one vial of monoclonal Anti-B reagent (clone numbers are mentioned on the corresponding certificate of analysis /release document and on the label). The reagents are prepared from culture supernatants from stable hybridoma cell lines as first described by Köhler and Milstein (Nature 1975). These monoclonal reagents contain murine IgM antibodies and have been specially selected and developed to provide a reliable alternative to polyclonal reagents. These reagents meet the requirements of the concerned standards and guidelines. Performance characteristics are mentioned in the certificate of quality control, which is supplied with every delivery. The principle of the test is the agglutination technique, which is based on antigen/antibody reaction. The reagents are used for examinations on diagnostic cards supplied with the set together with the sticks serving for mixing the reagents with blood samples.

Precautions

The diagnostic set is designed for professional in vitro use only. The set should be stored at 2 – 8°C. Leaking or damaged vials may not be used. Reagents (unopened or opened) should not be used beyond the expiration date, which is printed on the label of the vial. Sodium azide 0.1% (w/v) is used as preservative. The Anti-A and Anti-B reagents have been coloured for easy recognition. The reagents cannot be assumed to be free from infectious agents. Care must be taken in the use and disposal of each container and its contents. Turbidity may indicate microbial contamination. The disposal of waste after completion of the test should be performed according to local laboratory regulations.

Specimen collection and preparation

Blood samples should be taken aseptically.

Test procedure

Each test is done on one diagnostic card, the testing part of which is divided into two halves. The upper half is used to test the blood sample of the recipient - patient, the bottom half to test the blood sample of the donor-blood bag.

1. Instil one drop of the monoclonal Anti-A reagent in each blue circle on the card, and one drop of the monoclonal Anti-B reagent in each yellow circle.
2. Instil one drop of full blood of the recipient (patient) in the red circles in the upper half of the card, and one drop of the blood of the donor (blood bag) in the red circles in the bottom half of the card.

3. Mix the drops of monoclonal reagents with the drops of blood using the sticks supplied, each sample with a separate tip of the stick.

Interpretation

Read the results in 1 minute after stirring applying slight rocking motion with the diagnostic card. A positive reaction (i.e. agglutination) indicates the presence of the corresponding antigen on the erythrocytes. A negative reaction (i.e. no visible agglutination) indicates the absence of the corresponding antigen on the erythrocytes. Do not confuse a dry mix with agglutination!

The relevant blood group is determined from the following table:

Reagent reaction		Blood group
Anti-A	Anti-B	
+	-	A
-	+	B
+	+	AB
-	-	0

If the character of the reaction does not correspond with one of the above combinations, the cause of the different result must be identified before marking the blood group of the patient/donor. Blood is suitable for transfusion only if the results of the blood group determined in the upper half of the diagnostic card (patient) are identical to the results in the bottom half of the diagnostic card (donor).

Package

3 ml Anti-A reagent, 3 ml Anti-B reagent,
 30 diagnostic cards, 60 sticks.
 30 tests

Storage

Store at +2 – 8°C.

Limitations

Unexpected positive reactions may be due to: polyagglutination, pseudoagglutination, autoagglutination, mixed field reaction, the presence of Wharton's jelly together with umbilical cord cells. Unexpected negative or weak reactions may be due to weak antigens, mixed field reaction, decreased reactivity of the reagent. False positive or false negative results may be due to a test material contamination or to a deviation from the recommended technique.

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4. Daniels G.; Human Blood Groups. Blackwell Science Ltd. 1995.