











Sanquin Reagents B.V. Plesmanlaan 125 1066 CX Amsterdam The Netherlands	Phone: +31 20 5123599 Fax: +31 20 5123570 Reagents@sanquin.nl www.sanquin.org/reagents
--	---

LOT \_\_\_\_\_

Rh-Hr	Rh-Hr										Kell				Duffy		Kidd		Lewis		P	MNS				Luther		Xg			
	C	D	E	c	e	C <sup>w</sup>	f	*	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	P <sub>1</sub>	M	N	S	s	Lu <sup>a</sup>	Lu <sup>b</sup>	Xg <sup>a</sup>		
R <sub>1</sub> R <sub>1</sub>	0002410	+	+	0	0	+	0	/	/	+	+	0	+	0	+	+	+	0	0	+	+	+	+	0	+	0	+	0	+	0	Wr(a+)

Wra / 27012021

ID \_\_\_\_\_ LOT \_\_\_\_\_

Rh-Hr	Rh-Hr										Kell				Duffy		Kidd		Lewis		P	MNS				Luther		Xg		
	C	D	E	c	e	C <sup>w</sup>	f	*	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	P <sub>1</sub>	M	N	S	s	Lu <sup>a</sup>	Lu <sup>b</sup>	Xg <sup>a</sup>	
rr	0130896	0	0	0	+	+	0	/	/	0	+	0	+	/	+	0	+	+	+	w	+	0	+	0	0	+	0	+	+	Wr(a+)

Wra / 27012021

LOT \_\_\_\_\_

Rh-Hr	Rh-Hr										Kell				Duffy		Kidd		Lewis		P	MNS				Luther		Xg		
	C	D	E	c	e	C <sup>w</sup>	f	*	V	K	k	Kp <sup>a</sup>	Kp <sup>b</sup>	Js <sup>a</sup>	Js <sup>b</sup>	Fy <sup>a</sup>	Fy <sup>b</sup>	Jk <sup>a</sup>	Jk <sup>b</sup>	Le <sup>a</sup>	Le <sup>b</sup>	P <sub>1</sub>	M	N	S	s	Lu <sup>a</sup>	Lu <sup>b</sup>	Xg <sup>a</sup>	
rr	1283404	0	0	0	+	+	0	/	/	0	+	0	+	/	+	0	+	+	0	+	0	+	0	+	0	+	0	+	+	Wr(a+)

Wra / 27012021

+ = positive;положителен; pozitivní; positif; positif; θετική; positivo; positive; pozitív; positivo; pozitív; pozytywny; positivo; pozitiv; pozitivne; positif; pozitivno  
 o = negative; отрицателен; negativní; negativ; negativ; αρνητική; negativo; négative; negativ; negativo; negatív; negatywny; negativo; negativ; negativny; negativ; negativno  
 / = not tested; не е тестван; netestováno; ikke testet; nicht getestet; δεν έχει δοκιμαστεί; no analizado; non testée; nem tesztelt; non analizzato; nav pārbaudīts; nie testowane; não testado; netestat; netestované; inte testad; nije testirano  
 w = weak; слаб; slabý; svag; schwach; αδενής; débil; positive faible; gyenge; debole; vājš; slabý; fraco; slab; slabý; svag; slabo  
 s = strong; силно; silný; stærk; stark; ισχυρή; fuerte; positive forte; erősen; forte; spēcīga; silny; forte; puternic; silný; kraftig; jako  
 ID = donor number; брой на донорите; počet dárčů; donornummer; Donor-ID; αριθμός δότη; número de donante; numéro du donneur; donor száma; numero del donatore; donoru skaits; liczba dawców; número do dador; numărul de donator; počet darcov; donato