

Specifikace anti-D diagnostik

Clone	II	IIIa	IIIb	IIIc	IVa	IVb	Va	VI	VII	DFR	DBT	R ₀ ^{har}	HMi	Howe	Carr	Casar i	Lore	Sor	1711 844	
RUM-1	+	+	+	+	+	+	+	-	+	+	+	+	+	+	+	+	+	+	+	+
MS-201	+	+	+	+	+	+	+	-	+	+	+	+	+	+	+	+	+	+	+	+
TH-28	+	+	+	+	+	+	+	-	+	-	+	+	+	-	+	+	+	+	+	+
MS-26	-	+	+	+	-	-	+	+	+	+	-	-	+	+	+	+/-	+	+	+	+
MAD-2	+	NT	NT	NT	+	NT	(+)	-	+	-	-	-	-	NT	NT	NT	NT	NT	NT	NT