

Wr^a positive reagent red cells



CE 0344

IVD



Reagents

Plesmanlaan 125

1066 CX Amsterdam

The Netherlands

Phone: 31.20.512.3599

Fax: 31.20.512.3570

reagents@sanquin.nl

www.sanquinreagents.com

LOT _____ ⏳ _____

Rh-Hr	ID	Rh-Hr						Kell						Duffy		Kidd		Lewis		P					MNS					Luther.			Xg				
		C	D	E	c	e	C ^w	f	V*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a								
R ₁ r	M400409	+	+	0	+	+	0	/	0	+	+	0	+	0	+	+	+	0	+	0	0	+	+	+	+	0	+	+	+	+	0	+	+	+	+	Wr(a+)	Co(b+)

LOT _____ ⏳ _____

Rh-Hr	ID	Rh-Hr						Kell						Duffy		Kidd		Lewis		P					MNS					Luther.			Xg			
		C	D	E	c	e	C ^w	f	V*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a							
R ₁ r	M400233	+	+	0	+	+	0	/	0	+	+	0	+	0	+	+	+	0	0	+	+	+	+	+	+	0	+	+	+	+	+	0	+	+	+	Wr(a+)

LOT _____ ⏳ _____

Rh-Hr	ID	Rh-Hr						Kell						Duffy		Kidd		Lewis		P					MNS					Luther.			Xg			
		C	D	E	c	e	C ^w	f	V*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a							
R ₂ r	M500176	0	+	+	+	+	0	/	0	+	+	0	+	0	+	0	+	+	+	0	+	+	+	+	0	+	0	+	+	+	+	0	+	+	+	Wr(a+)

LOT _____ ⏳ _____

Rh-Hr	ID	Rh-Hr						Kell						Duffy		Kidd		Lewis		P					MNS					Luther.			Xg					
		C	D	E	c	e	C ^w	f	V*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a									
r'r	330982	+	0	0	+	+	0	/	/	0	+	0	+	0	+	+	+	+	0	0	+	0	+	0	+	0	+	0	+	+	+	+	0	+	0	+	Wr(a+)	Bg(a+)

Wr^a positive reagent red cells



CE 0344

IVD



Reagents

Plesmanlaan 125
1066 CX Amsterdam
The Netherlands

Phone: 31.20.512.3599

Fax: 31.20.512.3570

reagents@sanquin.nl

www.sanquinreagents.com

LOT _____



Rh-Hr	Rh-Hr									Kell				Duffy		Kidd		Lewis		P					MNS					Luther.			Xg
	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a				
rr	382612	0	0	0	+	+	0	/	/	+	+	0	+	0	+	+	+	+	0	0	+	0	+	+	0	+	0	+	+	+	+	Wr(a+)	

LOT _____



Rh-Hr	Rh-Hr									Kell				Duffy		Kidd		Lewis		P					MNS					Luther.			Xg
	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a				
r'r	M202322	+	0	0	+	+	0	/	0	0	+	0	+	0	+	+	+	+	0	+	0	+	+	0	+	+	0	+	+	+	Wr(a+)		

LOT _____



Rh-Hr	Rh-Hr									Kell				Duffy		Kidd		Lewis		P					MNS					Luther.			Xg
	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a				
rr	386169	0	0	0	+	+	0	/	/	+	+	0	+	0	+	+	+	+	0	0	+	+	+	+	+	+	0	+	+	+	Wr(a+)		

LOT _____



Rh-Hr	Rh-Hr									Kell				Duffy		Kidd		Lewis		P					MNS					Luther.			Xg
	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a				
R ₁ R ₁	002410	+	+	0	0	+	0	/	/	+	+	0	+	0	+	+	+	+	0	0	+	+	+	+	0	+	0	+	0	+	Wr(a+)		

Wr^a positive reagent red cells



CE 0344

IVD



Reagents
Plesmanlaan 125
1066 CX Amsterdam
The Netherlands

Phone: 31.20.512.3599
Fax: 31.20.512.3570
reagents@sanquin.nl
www.sanquinreagents.com

LOT _____



		Rh-Hr						Kell						Duffy		Kidd		Lewis		P				MNS				Luther.		Xg		
Rh-Hr	ID	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a		
rr	1372458	0	0	0	+	+	0	/	/	0	+	0	+	/	+	+	0	+	+	0	+	0	+	0	+	+	0	+	+	+	+	Wr(a+)

LOT _____



		Rh-Hr						Kell						Duffy		Kidd		Lewis		P				MNS				Luther.		Xg	
Rh-Hr	ID	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	
rr	337475	0	0	0	+	+	0	/	/	+	+	0	+	/	+	+	0	0	+	+	0	0	0	+	+	0	+	0	+	+	Wr(a+)

LOT _____



		Rh-Hr						Kell						Duffy		Kidd		Lewis		P				MNS				Luther.		Xg	
Rh-Hr	ID	C	D	E	c	e	C ^w	f	V	*	K	k	Kp ^a	Kp ^b	Js ^a *	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	
rr	130896	0	0	0	+	+	0	/	/	0	+	0	+	/	+	0	+	+	+	W	+	0	+	0	0	0	+	0	+	+	Wr(a+)

Wra / 010312

+ = positive; positif; positive; positivo; positivo; positivo; θετική; positiv; positiv; pozytywny; pozitivni; pozitivne; pozitiv

0 = negative; negativ; négative; negativo; negativo; negativo; αρνητική; negativ; negativ; negatywny; negativni; negativny; negativ

/ = not tested; nicht getestet; non testé; non analizzato; no analizado; não testado; δεν έχει δοκιμαστεί; inte testad; ikke testet; nie testowane; netestováno; netestované; nem tesztelt

W = weak; schwach; positive faible; debole; débil; fraco; ασθενής; svag; svag; slaby; slabý; slabý; gyenge

S = strong; stark; positive forte; forte; fuerte; forte; ισχυρή; kraftig; stærk; silny; silný; silný; erős

ID = donornumber; Donor-ID; numéro du donneur; numero del donatore; número de donante; número do dador; αριθμός δότη; donatornummer; donornummer; liczba dawców; počet darců; počet darcov; donor száma